# Decode-seq3

TuLab IGDB CAS, 2024-01

# **Overview**

- Decode-seq1 protocol: http://tulab.genetics.ac.cn/media/filer\_public/76/35/763595a7-0ea6-46ed-9220-9e38717ec384/decode-seq\_protocol.pdf
- prime-seq protocol: https://www.protocols.io/view/prime-seq-81wgb1pw3vpk/v2
- Decode-seq3 是在 Decode-seq1 的基础上做出适当修改,主要参考了 prime-seq 在 RNA 提取、 反转录酶使用、建库等方面的改进。三种方法的主要区别如下图所示:



• Decode-seq3 文库构建流程如下图所示:

5'-RNA:NB(A)<sub>30</sub>- 3'

3'-CCC:cDNA:NV(T)<sub>30</sub> GCTGA CGCAGCACATCCCTTTCTCACA-5' (5bc-RT Primer)

Template switch (5bc-TSO; 5bc-RTb, b 表示加了 Biotin) 5'-ACACTCTTTCCCTACACGACGCTCTTCCGATCT[BC<sub>6</sub>](N)<sub>17</sub>GGG:RNA:NB(A)<sub>30</sub>- 3' (PE Read 1 Sequencing Primer) 3'-CCC: cDNA: NV (T) 30 GCTGA CGCAGCACATCCCTTTCTCACA- 5' 3' -<u>TGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA</u>[BC<sub>6</sub>] (N) 17CCC:cDNA:NV(T) 30 GCTGA <u>CGCAGCACATCCCTTTCTCACA</u>- 5' Preamp(single-Pb, 同 pm-Preamp) 3'-CGCAGCACATCCCTTTCTCACA-5' (single Pb) 5' -ACACTCTTTCCCTACACGACGCTCTTCCCGATCT [BC6] (N) 17GGG: cDNA:NB (A) 30 CGACT GCGTCGTGTAGGGAAAGAGTGT- 3' 3' -<u>TGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA</u>[BC<sub>6</sub>] (N) 17CCC:cDNA:NV (T) 30 GCTGA <u>CGCAGCACATCCCTTTCTCACA</u>- 5' 5'-<u>ACACTCTTTCCCTACACGACGC</u>-3' (single Pb) 使用 NEB 试剂盒(最终采用) **Fragmentation** P-ACACTCTTTCCCTACACGACGCTCTTCCGATCT[BC6](N)17GGG:5'Frag-A A-TGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA[BC<sub>6</sub>](N)<sub>17</sub>CCC:5'Frag-P Ligation (pm-Adap) GTCTCGTGGGGCTCGGAGATGTGTATAAGAGACAGT P-ACACTCTTTCCCTACAGGACGCTCTTCCGATCT[BC6](N)17GGG:5'Frag-A P- CTGTCTCTTATACACATCT TCTACACATATTCTCTGTC-P A-TGTGAGAAAGGGATGTGCCGGGGAAGGCTAGA[BC<sub>6</sub>](N)<sub>17</sub>CCC:5'Frag-P TGACAGAGAATATGTGTAGAGGCTCGGGTGCTCTG GTCTCGTGGGGCTCGGAGATGTGTATAAGAGACAGTACACTCTTTCCCTACACGACGCTCTTCCGATCT[BC<sub>6</sub>](N)<sub>17</sub>GGG:5'Frag-ACTGTCTCTTATACACATCT TCTACACATATTCTCTGTCATGTGAGAAAGGGATGTGCTGCGAGAAGGCTAGA[BC6](N)17CCC:5'Frag-TGACAGAGAATATGTGTAGAGGCTCGGGTGCTCTG <mark>PCR (5' enrich)</mark> 理论上不丢失任何 5'reads (P5R1, 同 pm-3enrich; i7 index primer) AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGTACACTCTTTCCCTACACGACGCTCTTCCGATCT[BC6](N)17GGG:5'Frag-ACTGTCTCTTATACACATCT TCTACACATATTCTCTGTCATGTGAGAAAGGGATGTGCCGGGGAGGCTGGA[BC6](N)17CCC:5'Frag-TGACAGAGAATATGTGTAGAGGCTCGGGTGCTCTG (Nextera i7 Index) 3'-GGCTCGGGTGCTCTG[i7]TAGAGCATACGGCAGAAGACGAAC-5' Library

AATGATACGGCGACCACCGAGATCTACAC<mark>TCTTTCCCTACACGACGC</mark>TCTTCCGATCT[BC6](N)<sub>17</sub>GG:5'Frag-ACTGTCTCTATACACATCTCCGAGCCCACGAGAC[i7]ATCTCGTATGCCGTCTTCTGCTTG TTACTATGCCGCTGGTGGCTCTAGATGTG<mark>AGAAAGGGATGTGCTGCGAGAAGGCTAGA[</mark>BC6](N)<sub>17</sub>CC:5'Frag-TGACAGAGAATATGTGTAGAGGCTCGGGTGCTCTC[i7]TAGAGCATACGGCAGAAGACGAAC-5'

## Preparation

### 试剂购买及配制

#### 试剂购买

name	cat
Sera-Mag Speed Beads	Ge Healthcare cat. 65152105050250

name	cat	
PEG 8000, Poly(ethylene glycol)	Sigma Aldrich cat. 89510	
EDTA (0.5M, PH8.0)	BioRoYee cat. N-R1007,	
Igepal(10%) 非离子型表面活性剂	Sigma Aldrich cat. 18896	
TE Buffer, PH8.0	Life Technologies cat. AM9849	
NaCl (5 M), 无 RNase	Invitrogen cat. AM9760G	
Tris-HCI (1M, pH 8.0)	BioRoYee cat. NR0073	
Maxima H Minus Reverse Transcriptase	Thermo Scientific cat. EP0753	
dNTPs (10 mM)	Thermo Scientific cat. R0192	
Exonuclease I (E.coli)	NEB cat.M0293V	
KAPA HiFi 2x ReadyMix	KAPA cat. kk2602	
NEBNext Ultra II FS DNA Library Prep Kit for Illumina	NEB cat. E7805S	
Buffer EB	Qiagen cat. 19086	
SPRI select beads	Beckman Coulter, cat. B23317	

• 📓 Ref: Decode-seq1 protocol 和 prime-seq protocol (见 Overview 部分)

#### 试剂配制

- **PEG Solution** (22%)
  - Prepare PEG Solution (22%) by adding all ingredients to a 50 mL falcon tube. Do not add the total amount of water until after PEG is completely solubilized.
  - Incubate at 40 °C and vortex regularly until PEG is completely dissolved
  - Notes: 最终 49ml PEG Solution 中加入了1ml beads,相当于beads浓度为 2%,应该足够 吸附DNA。

Reagent	Amount
PEG 8000	11 g
NaCI (5M)	10 mL
Tris-HCI (1M, pH 8.0)	500 μL

Reagent	Amount	
EDTA (0.5M)	100 μL	
IGEPAL (10% solution)	50µL	
UltraPure Water	up to 49 mL	
Total	49mL	

- Cleanup Beads (in 22% PEG Solution)
  - i. Resuspend **Sera-Mag Speed Beads** carefully and pipette **1 ml** of bead suspension into a 1.5 mL tube
  - ii. Place on magnet stand and remove storage buffer
  - iii. Add 1 ml of TE Buffer and resuspend beads
  - iv. Place on magnet stand and remove supernatant
  - v. Repeat wash step one more time
  - vi. Add 900  $\mu I$  TE Buffer and resuspend beads

vii. Add the washed Sera-Mag Speed Beads to the 49mL PEG Solution (22%) and mix well

- Annealing buffer (10mM Tris pH8, 50mM NaCl, 1 mM EDTA)
  - 990µL TE + 10µL 5M NaCl

### **Primer sequence**

• primers: thaw on ice

Primer name	Primer sequence
5bc- RTb	Biotin-ACACTCTTTCCCTACACGACGC AGTCG T30 VN
5bc- TSO	ACACTCTTTCCCTACACGACGCTCTTCCGATCT BC6 N17 /rG/rG/rG/
single- Pb(同 pm- Preamp)	Biotin-ACACTCTTTCCCTACACGACGC

Primer name	Primer sequence
pm- Adapter- Se	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGT
pm- Adapter- Anse	/5Phos/CTGTCTTATACACATCT
prime- seq Adapter	pm-Adapter-Se 和 pm-Adapter-Anse 的退火产物
P5R1(同 pm- 3enrich)	AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT
i7 index primer	CAAGCAGAAGACGGCATACGAGAT[i7]GTCTCGTGGGCTCGG

- 退火步骤:
  - 用 Annealing buffer 稀释引物(pm-Adapter-Se 和 pm-Adapter-Anse)至100µM,分别取 30µL,等量混合;
  - 。在PCR仪上运行退火程序;

		Cycles	Temperature	Time
Advanced Protocol	Step 1:	1	95°C	5 min
(example in which the	Step 2:	40*	95°C (-1°C/cycle)	1 min
oligonucleotide pair	Step 3:	1	55°C	30 min
has a T <sub>m</sub> of 55°C)*	Step 4:	20*	55°C (-1°C/cycle)	1 min
	Step 5:		4°C	HOLD

\*The number of cycles in Step 2 and 4 depends on the  $T_m$  of the oligonucleotides to be annealed.

。 退火产物 prime-seq Adapter 此时浓度为 50μM, 稀释至 1.5μM 后使用。

# **RNA** extraction

- RNA 提取可自行选择方法(Trizol、Qiagen柱子、beads)
- Beads 提取 RNA 主要包括下列步骤(参考 prime-seq protocol):
  - Sample Collection
  - Proteinase K Digest

- Bead Clean Up
- DNase I Digest

# **Reverse Transcription (Decode3)**

• Prepare RT Mix:

Reagent	Well	Plate
Maxima H Minus RT	0.15 µL	16.5 μL
Maxima RT Buffer (5x)	2 µL	220 µL
dNTPs (10 mM)	1 μL	44*2.5 μL
5bc-RTb primer(10uM)	1 μL	
Total	4.15 µL	

- add RNA samples per well (4-40ng RNA), final volumn 5.75 µL using RF water.
- Add 0.1 μL **5bc-TSO** (100 μM) per well.
  - 实验室有 30 种 5bc-TSO,其中 5bc-TSO30,5bc-TSO35,5bc-TSO38 测试效果不好,分别用 5bc-TSO5,5bc-TSO6,5bc-TSO7 代替。

```
5bc-TS011, 5bc-TS012, 5bc-TS013, 5bc-TS014, 5bc-TS015,
5bc-TS016, 5bc-TS017, 5bc-TS018, 5bc-TS019, 5bc-TS020,
5bc-TS021, 5bc-TS022, 5bc-TS023, 5bc-TS024, 5bc-TS025,
5bc-TS026, 5bc-TS027, 5bc-TS028, 5bc-TS029, 5bc-TS030(TS05),
5bc-TS031, 5bc-TS032, 5bc-TS033, 5bc-TS034, 5bc-TS035(TS06),
5bc-TS036, 5bc-TS037, 5bc-TS038(TS07), 5bc-TS039, 5bc-TS040
```

Incubate for 1:30 at 42°C

# cDNA pooling & beads purification (Decode3/pm)

- Pool all wells cDNA into a 1.5/2 mL tube
- Add 10 µl of Cleanup Beads (22% PEG) for each sample for a 1:1 ratio

- 📓 Notes: The EDTA in the Cleanup Beads will inactivate the RT)
- ■ Notes: Final PEG 11%, 1.8x Ampure beads 为 PEG 18%\*1.8/2.8 = 11.57%, left side size 约为 150bp(参考 beads 测试电泳图)
- Incubate for 5min at RT to allow binding of the cDNA onto beads
- Place the tube on the magnet stand until clear (~3 min) and discard supernatant
- Wash with **1 mL** of **80% EtOH** while the tube is on the magnet. Discard the supernatant.
  - Notes: Volume of EtOH should be adjusted depending on the number of samples. More samples will require more EtOH to cover the beads completely.
  - I When washing the sample with 80% ethanol, keep the sample tube on the magnetic frame and do not stir the magnetic beads.
  - I The volume of 80% ethanol should be larger than the original RNA and beads mix
- Repeat wash step once more
- Air dry beads for 5min (regularly check the beads and avoid over-drying)
- Elute the beads in 17 µL of UltraPure Water
- Incubate for **5min** at RT and transfer to a new PCR tube or plate.

# Exonuclease I Treatment (Decode3/pm)

• Add 2 µI of Exol Buffer (10x) and 1 µI of Exonuclease I. Incubate as follows:

Step	Temperature	Time
Incubation	37°C	20 min
Heat Inactivation	80°C	10 min
Storage	4°C	hold

- Mix each sample (20 µL per well) with 16 µl of Cleanup Beads (22% PEG) for a 1:0.8 ration
  - INOTES: final PEG 9.778%, 1.2x anpure beads 为 PEG 18%\*1.2/2.2 = 9.82%, left side size 约为 200-250bp
- Incubate for 5min at Room temperature to allow binding of the cDNA onto beads
- Place the tube on the magnet stand until clear (~3 min) and discard supernatant
- Wash with 50  $\mu l$  of 80% EtOH while the tube is on the magnet. Discard the supernatant
- Repeat wash step once more
- Air dry beads for 5min
- Elute the beads in 20 µl of UltraPure Water
- Incubate for 5min at RT and transfer to a new PCR tube or plate

# Full length cDNA Amplification (Decode3/pm)

• Prepare Pre Amplification Mix

Reagent	1x
KAPA HiFi 2x	25 μL
single-Pb (10uM)	3 µL
NF Water	2 µL
Total	30 μL

- Add **30 µI** Pre Amplification Mix to sample
- Incubate the Pre Amplification PCR as follows:

Step	Temperature	Time	Cycles
Initial Denaturation	98°C	3 min	1 cycle
Denaturation	98°C	15 sec	10 cycles*
Annealing	65°C	30 sec	
Elongation	72°C	4 min	
Final Elongation	72°C	10 min	1 cycle
Storage	4°C	hold	

• Adjust the number of cycles based on input (sample number, cell number, or concentration). As a general guide we recommend:

Total RNA Input	Cycles
10 ng	16
50 ng	14
100 ng	12
500 ng	10

# cDNA Bead Purification (Decode3/pm)

- Mix sample with 40µI Cleanup Beads (22% PEG) for a ratio of 1:0.8
- Incubate for 5min at 20°C (Room Temp)
- Place the tube on the magnet stand until clear (~3 min) and discard supernatant
- Wash with 100  $\mu I$  of 80% EtOH while the tube is on the magnet. Discard the supernatant
- Repeat wash step once more
- Air dry beads for 5min
- Elute cDNA in 10 µl UltraPure Water
- Incubate for 5min at RT and transfer to a new PCR tube or plate

Stopping Point. Samples can be safely stored at -20°C and protocol can be continued at a later date.

- Quantify the cDNA using Qubit. Use **1 µI** of clean cDNA for quantification.
- Quality check the cDNA using the Agilent 2100 Bioanalyzer with High Sensitivity DNA Analysis Kits. 或者跑胶看一眼。
- 产物够多次建库的量,稀释至 4-8 ng/µL 用于后续 NEB protocol 或者 0.2 ng/µL 用于 Nextera protocol

### Library Preparation Using NEB (Decode3/pm)

### **NEB Fragmentation**

• Prepare Fragmentation Mix

Reagent	1x
Ultra II FS Reaction Buffer	1.4 µL
Ultra II FS Enzyme Mix	0.4 µL
cDNA (4-8 ng/µL)	2.5 μL
TE	1.7 µL

Reagent	1x
Total	6 µL

- Notes: Ensure that the Ultra II FS Reaction Buffer is completely thawed. If a precipitate is seen in the buffer, pipette up and down several times to break it up, and quickly vortex to mix.
   Place on ice until use.
- Notes: Vortex the Ultra II FS Enzyme Mix for 5-8 seconds prior to use for optimal performance.
- Vortex the Fragmentation Mix for 5s and immediately proceed to next step.
- Incubate the Fragmentation reaction as follows:

Step	Temperature	Time
Pre-Cool	4 C	×
Fragmentation	37 C	5 min
A Tailing and Phosphorylation	65 C	30 min
Storage	4 C	∞

• 📓 Notes: Set heated lid to 75° C.

### **NEB Adapter Ligation**

• Prepare Adapter Ligation Mix

Reagent	1x
NEBNext Ultra II Ligation Master Mix	6 µL
NEBNext Ligation Enhancer	0.2 µL
prime-seq Adapter (1.5 µM)	0.5 µL
Total	6.7 μL

- Add 6.7 µI Adapter Ligation Mix to each replicate
- Incubate for 15min at 20°C
  - Notes: Turn off heated lid
- Add 37.3 µI Buffer EB to Samples

- Buffer EB QIAGEN: The composition of Buffer EB is: 10 mM Tris-Cl, pH 8.5。暂用 pH 8.0 Tris-Cl 代替
- Mix Index PCR with 26 µl SPRI select beads (ratio of 1:0.52)
  - $\circ$   $\$ Notes: 26/50=0.72
  - Notes: We use SPRI Select Beads here instead of our home made 22% Cleanup beads for their guaranteed QCed size selection properties.
- Incubate for 5min at 20°C (Room Temp)
- Place the plate on the magnet stand until clear and transfer 76 µl supernatant to clean well.
- Mix supernatant with 10 µl SPRI select beads (ratio of 1:0.72)
  - Solution Notes: (26 + 10/50) = 0.72
- Incubate for 00:05:00 at 20°C (Room Temp)
- Place the plate on the magnet stand until clear and discard supernatant.
- Wash with 150 µl of 80% EtOH while the plate is on the magnet. Discard the supernatant
- Repeat wash step once more
- Air dry beads for 5min
- Elute samples in 10.5 µl 0.1X TE (dilute 1X TE Buffer 1:10 in water) for 5min

### **NEB Library PCR**

- Transfer samples to clean wells
- Add 1µl of i7 Index Primer (5 uM) to each well
  - 📓 Notes: This is the unique index that will be used for demultiplexing libraries.
- Prepare Library PCR Mix

Reagent	1x	5x
NEBNext Ultra II Q5 Master Mix	12.5 µL	62.5 μL
P5R1/pm-3enrich (5uM) 🚖	1 µL	5 µL
Total	13.5 µL	67.5 μL

- Notes: Although scaled down, there will not be sufficient Q5 Master Mix (M0544S 1.25ml /M0544L 6.25ml) in the kit. This item will have to be ordered separately.
- SAQ: What is the difference between NEBNext Ultra II Q5 Master Mix and NEBNext High-Fidelity 2X PCR Master Mix (M0541)? Although both master mixes contain Q5 High-Fidelity DNA Polymerase, the master mix composition is different. NEBNext Ultra II Q5 Master Mix has been specifically optimized for to faithfully amplify NGS libraries with the least amount of amplification bias. Additionally, amplification of NGS libraries in the presence or absence of magnetic beads is more robust, and it is a hot start, making it ideal for workflows involving automation.

- Add **13.5µl** of Library PCR Mix to each well
- Incubate the Library PCR reaction as follows:

Step	Temperature	Time	Cycles
Initial Denaturation	98 C	30 sec	1 cycle
Denaturation	98 C	10 sec	10 cycles*
Annealing/Elongation	65 C	1 min 15 sec	
Final Elongation	65 C	5 min	1 cycle
Storage	4 C	hold	

• Adjust the number of cycles based on cDNA input. As a general guide we recommend.

cDNA Input	Cycles
20 ng	10
10 ng	11
5 ng	12

### Library Double Size Selection (Decode3/pm)

- Add 25 µl Buffer EB to Index PCR (补齐PCR体系体积至50µL)
- Mix Index PCR with 26 µl SPRI select beads (ratio of 1:0.52)
- Incubate for 00:05:00 at 20°C (Room Temp)
- Place the plate on the magnet stand until clear and transfer **76 µI** supernatant to clean well.
- Mix supernatant with 10 µl SPRI select beads (ratio of 1:0.72)
- Incubate for 5min at 20°C (Room Temp).
- Place the plate on the magnet stand until clear and discard supernatant.
- Wash with 150 µl of 80% EtOH while the plate is on the magnet. Discard the supernatant.
- Repeat wash step once more.
- Air dry beads for 5min.
- Elute in **15 µI** UltraPure Water.
- Incubate for 5min and then place on magnet until clear. Transfer eluted library to new well.